

ECO TRAVEL

YOUR TICKET TO EFFORTLESS MOBILITY

- FOLDABLE HANDLEBARS
- DUAL BREAKING SYSTEM
- 5 COLOURS
- SMART CASE AVAILABLE*



* check available options and models in the table



Rehasense's quality management system complies with both the ISO 9001 and ISO 13485 standards, as confirmed by TUV NORD. Rehasense's medical devices are CE-marked and meet the requirements of EU Regulation 2017/745 (MDR). Rehasense's medical devices are UKCA-marked and meet the requirements of MHRA regulations.



Errors and omissions excepted

RMFTALLENO

ECO TRAVEL

THE SMART WAY TO TRAVEL



REHASENSE

NEW
PRODUCT



www.rehasense.com

ECO TRAVEL

COMPACT & TRAVEL-FRIENDLY

Introducing the Rehasense ECO TRAVEL power add-on for manual wheelchairs - a game-changer for those seeking enhanced mobility on the go. Designed with your mobility in mind, the ECO TRAVEL is the lightest power add-on in our lineup, compact and easy to take with you wherever you go. Whether you are at home, on the road, or exploring new places, it has been built for convenience and independence. The ECO TRAVEL effortlessly transforms your manual wheelchair into a powerful, easy-to-manoeuvre companion, perfect for both short trips and extended journeys. You can choose between the two versions ECO TRAVEL STANDARD (Bx) or ECO TRAVEL ADVANCED (Lx).

EASY SETUP & PERFORMANCE

- **Quick & easy installation:** the modular mounting system requires no pre-mounted brackets, allowing for hassle-free setup
- **Customizable comfort:** user-friendly settings, including height and length adjustments, ensure optimal sitting and driving comfort
- **Unmatched battery range:** longest battery range on the market in this category - up to 45 km with standard battery
- **Stable & secure:** the 2-clamp attachment on the front frame provides optimal driving stability and protection

COMPATIBILITY

- Can be used with rigid and folding manual wheelchairs - and also models with removable leg rests
- Can be mounted on almost all wheelchair models on the market*
- Covers mounting frame widths between 220-480 mm


DOCKING & UNDOCKING

Swift Docking and Undocking: Experience the fastest setup on the market - quick and effortless, with only seconds needed to be ready to drive

SAFETY YOU CAN TRUST

- Dual braking system (drum and electronic brakes)
- The self-centering system makes driving safe and comfortable

*We do not recommend using ECO TRAVEL with carbon fibre wheelchairs. Rehasense takes no responsibility for any damage.



SPECIAL FEATURES



Dual braking system (drum and electronic brakes)



TECHNICAL DATA

	ECO TRAVEL Bx (STANDARD)	ECO TRAVEL Lx (ADVANCED)
Total weight	13.08 kg	14.16 kg
Transport weight	8 kg	9.04 kg
Max. user weight	110 kg	
Wheel size	8"/200 mm	
Engine power	250 W	350 W
Max. engine torque	16 Nm	19 Nm
Speed modes	2 gears	3 gears
Max speed	up to 15 km/h	up to 20 km/h
Climbing power	up to 10% (6°)	
Standard battery	Li-on, 48 V/ 5,8 Ah/ 278 Wh	
Driving range	up to 45 km	
Display	1.34" LCD black/white	1.81" LCD, color
· Setting options	None	Yes
· Language	English only	Multilingual
· Cruise control	No	Yes
· Motor brake (CLC)	No	Yes
· WALK function	No	Yes
Colours/Personalisation	Black	Green, Orange, Blue, Red, Black
Handlebar	Fixed	Foldable
Handlebar stem	100 mm, fixed length	Length adjustable 180-215 mm and easily foldable (without tools)
Support legs type	Simple feet	Incl. wheels
Travel suitcase	Available as an accessory	Included

We reserve the right to make changes.



Introducing the ECO TRAVEL

Smart Case! With a custom-fit interior, shock-absorbing padding, and lightweight, canvas material, it makes moving your ECO TRAVEL safe, simple, and secure - whether you are traveling near or far.

Height: 77 cm
Width: 51 cm
Depth: 36 cm
Travel bag incl. EPE: 8.6 kg
Charger and tools: 1 kg

Custom-fit interior

