

# ICON 40 X

THE ADAPTIVE WHEELCHAIR



REHASENSE



# ICON 40 X

## THE ADAPTIVE WHEELCHAIR



**COLOURS**  
Champagne, Silver



The ICON 40 is the most adjustable model within the BASE LINE range. It has a seat depth adjustable frame with multiple wheels, axle, cross bar and front caster set up options. This ranges from active self-propelling drivers to rehabilitation purposes, where a high degree of control, support and stability is required.

### FEATURES AT A GLANCE

- Solid aluminium wheelchair
- Adjustable depth frame
- Swing-away, fixed height, detachable arm rest
- Adjustable leg rests via bolts: swing in and out, removable
- Height-adjustable leg rests via bolts or clamps
- Flip-up footplate
- Nylon back rest and seat for superb comfort and fit
- Wide choice of frame sizes available to ensure every user has the correct posture and best driving experience
- Quick-release axle
- Multiple options for anti-tip: slide up version, swing away version or swing up version
- 24" PU rear wheel and 8" PU front wheel

### TECHNICAL DETAILS

Weight from	Transport weight from	SWL (Safe Weight Limit)	Seat widths	Seat depth	Total width	Total length	Total width folded	Total height	Backrest height	Armrest height
15.4 kg	8,0 kg	150 kg	40 cm (standard), further options: 30, 35, 37.5, 42.5, 45, 47.5, 50, 53, 57, 61	30-45 S* / 40-55 L* cm	Seat width + 21 cm	101 cm	34 cm	92,5 cm	42 cm	23-28 cm
		190 kg								

\* S short frame, L long frame

### SPECIAL FEATURES



Detachable and height-adjustable armrests with side guards that swing back for transfers



Legrests "swing in, swing out, swing under" as well as removable. Flip-up, angle-adjustable foot plates with heel band



Solid puncture proof rear wheel with push rim. Easy to disconnect with the quick-release axle




Adjustable front fork (angle and height) to align perpendicular to the ground for optimal driving




Slide up height adjustable anti-tip combined with step assist for easier tipping driving




Standard push operation hand brakes



Multiple rear wheel positions for variation in seat depth and seat height



Height adjustable push handles via back cane



Your Rehasense Dealer

VI-2022

Rehasense's quality management complies with both the ISO 9001 and ISO 13485 standards, as confirmed by TÜV NORD. Rehasense's medical devices are CE-marked and meet the requirements of EU Regulation 2017/745 (MDR).

